

**OXIDIZED FUNGAL ANTIGENS AND METHODS OF MAKING AND USING
THEREOF**

Technical Field

The present invention relates to oxidized fungal antigens and methods of making
5 and using thereof. More particularly, the present invention provides a method for
producing an oxidized fungal antigen in culture filtrate. The present invention also
provides for the produced oxidized fungal antigens. Devices comprising such oxidized
fungal antigens, methods for testing for fungal antibodies using the oxidized fungal
antigens and methods for producing anti-fungal antibodies using oxidized fungal antigens
10 are further provided. Antigen detection devices comprising anti-fungal antibodies raised
against oxidized fungal antigens produced by the present methods are further provided.

Background Art

Fungal infections are common in the United States, and are often associated with
15 significant morbidity. For example, *Candida* species are the fourth most common
organisms isolated from US hospital patients (Wenzel and Pfaller, *Infect. Control Hosp.*
Epidemiol., 12:523-524 (1991)). *Histoplasma capsulatum* and *Coccidioides immitis*
causes, respectively, 250,000 and 100,000 new infections annually, and 10-40% of
infected persons become symptomatic (Bullock, *Histoplasma capsulatum*, In: *Principles*
20 and *Practice of Infectious Diseases*, GL Mandell et al. (eds.), Churchill Livingstone, New
York, 1995, pp. 2340-2353 and 2365-2375; and Stevens, *N. Engl. J. Med.*, 332:1077-
1082 (1995)). Individuals with symptomatic *C. immitis* infection miss an average of 35
days of school or work, and incur a total of \$24 million in medical expenses annually.

Many fungi cause respiratory infections that are indistinguishable on the basis of
25 symptoms. However, it is important to accurately determine which fungus is responsible
for infection due to differences in treatment regimens and incidences of complications
(Bullock, *Histoplasma capsulatum*, In: *Principles and Practice of Infectious Diseases*, GL
Mandell et al. (eds.), Churchill Livingstone, New York, 1995, pp. 2340-2353 and 2365-
2375; and Galgiani, *Ann. Intern. Med.*, 130:293-300 (1999)). Because culture and direct
30 antigen detection are problematic, antibody detection plays an important role in the

diagnosis and identification of fungal infections. The methods most commonly used to detect fungal antibodies include complement fixation (CF) and immunodiffusion (ID) (Zancope-Oliveira et al., *Clin. Diagn. Lab. Immunol.*, 1:90-93 (1994)). Both methods, however, have inherent drawbacks. ID detects antibodies to species-specific protein
5 moieties, and is thus sensitive and highly specific; however, the assay takes 48 hours to perform, and appropriate interpretation requires highly skilled personnel. Although CF assays are highly sensitive, their performance is complex and labor-intensive, and they exhibit low specificity due to cross-reactive antibodies recognizing carbohydrate moieties common to several fungi (Zancope-Oliveira et al., *Clin. Diagn. Lab. Immunol.*, 1:90-93
10 (1994); and Yang et al., *Clin. Diagn. Lab. Immunol.*, 4:19-22 (1997)). Enzyme-linked immunosorbent assays (ELISA) for the detection of fungal antibodies have also been described, but these assays exhibit low specificity due to the same cross-reactive
15 antibodies at issue in CF assays (Zartarian et al., *Am. J. Clin. Pathol.*, 207:148-153 (1997); and Kaufman et al., *J. Clin. Microbiol.*, 33:618-619 (1995)).

During the last decade, investigators using a western blot assay system found that periodate oxidation of fungal antigen preparations inactivates the cross-reactive carbohydrate moieties but does not disturb the structural integrity of the protein moieties (Zancope-Oliveira et al., *Clin. Diagn. Lab. Immunol.*, 1:90-93 (1994); Zancope-Oliveira et al., *Clin. Diagn. Lab. Immunol.*, 1:563-568 (1994); and Pizzini et al., *Clin. Diagn. Lab. Immunol.*, 6:20-23 (1999)). Thus, utilization of oxidized fungal antigens markedly increased the specificity of the western blot assay for fungal antibodies, without affecting sensitivity. However, the western blot assay described therein is not practical for most medical laboratories; many complex steps are required, including antigen preparation, gel electrophoresis, electrotransfer to nitrocellulose membranes, assay performance, and blot
20 interpretation. The antigen was a culture filtrate, 0.45 μ m filtered, 20x concentrated by ultrafiltration, dialyzed against PBS, purified by ion exchange chromatography, 10X concentrated by ultrafiltration, and finally dialyzed against PBS. Oxidizing the antigen also included several labor intensive and potentially contaminating steps: 4mg/mL antigen was exposed to periodate for 18 hours at 4°C, treated with glycerol for 15
25 minutes, treated with borohydride for 2 hours, and finally dialyzed against water. Manufacturing the devices also included methods requiring sophisticated equipment and
30

skilled technicians: denaturing the oxidized antigen with heat and mercaptoethanol, electrophoresing with a 4% stacking gel and 7.5% resolving gel, transferring the proteins to nitrocellulose paper, blocking with buffered milk and tween, drying the paper, and cutting the paper into strips, etc.

5 Although attempts to utilize periodate-treated fungal antigens for specific detection of fungal antibodies in an ELISA system have been reported, these attempts were only marginally successful (Fisher et al., *Mycoses*, 40:83-90 (1997)). Binding of *H. capsulatum* antibodies to oxidized *Blastomyces dermatitidis* antigens was reduced only by about 20% when compared to binding to non-oxidized *B. dermatitidis* antigens.

10 This limited reduction in cross-reactivity may have been due to the use of highly purified (and thus extensively manipulated) fungal antigens, periodate oxidation after the antigens had adsorbed to the plastic surface, or the use of an insufficient concentration of periodate.

15 Therefore, there exists a need in the art for oxidized fungal antigens with more reduced cross-reactivity. The present invention addresses this and other related needs.

Disclosure of the Invention

20 In one aspect, the present invention is directed to a method for producing a fungal antigen suitable for testing for an antibody to a fungus, the method comprising: a) providing a fungal antigen as a culture filtrate, the fungal antigen having not been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate; and b) contacting the fungal antigen with an oxidizing agent to produce an oxidized fungal antigen suitable for testing for an antibody to the fungus. Devices comprising an oxidized fungal antigen, or a plurality of oxidized fungal antigens, produced by the above 25 method, for testing for an anti-fungal antibody, or a plurality of anti-fungal antibodies, are also provided.

30 In another aspect, the present invention is directed to a method of testing for an antibody to a fungus in a sample, the method comprising: a) producing a fungal antigen suitable for testing for an antibody to a fungus, comprising providing a fungal antigen as a culture filtrate, the fungal antigen having not been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate, and contacting the fungal

antigen with an oxidizing agent to produce an oxidized fungal antigen suitable for testing for an antibody to the fungus; b) contacting a sample suspected of containing an antibody to a fungus with the oxidized fungal antigen produced in step a) under suitable conditions to allow binding of the antibody, if present in the sample, to the oxidized fungal antigen; 5 and c) assessing binding between the antibody and the oxidized fungal antigen to determine the presence and/or amount of the antibody in the sample.

In still another aspect, the present invention is directed to a method for producing an antibody to a fungal antigen, the method comprising: a) producing a fungal antigen comprising providing a fungal antigen as a culture filtrate, the fungal antigen having not 10 been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate, and contacting the fungal antigen with an oxidizing agent to produce an oxidized fungal antigen; b) delivering, to a vertebrate, the oxidized fungal antigen, in an amount sufficient to induce detectable production of an antibody to the antigen; and c) recovering the antibody from the vertebrate.

15 In yet another aspect, the present invention is directed to a method for producing a monoclonal antibody to a fungal antigen, the method comprising: a) producing a fungal antigen comprising providing a fungal antigen as a culture filtrate, the fungal antigen having not been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate, and contacting the fungal antigen with an oxidizing agent to produce an oxidized fungal antigen; b) delivering, to a vertebrate, the oxidized fungal antigen, in an amount sufficient to induce detectable production of an antibody to the antigen; c) removing at least a portion of antibody-producing cells from the vertebrate; d) immortalizing the removed antibody-producing cells; e) propagating the immortalized antibody-producing cells; and f) harvesting monoclonal antibody produced by the 20 immortalized antibody-producing cells. Methods of testing for a fungal antigen in a sample using the produced antibodies are also provided.

Brief Description of the Drawings

30 Figure 1 illustrates a *Histoplasma capsulatum* IgG ELISA test using an oxidized and a non-oxidized antigen.

Figure 2 illustrates a *Coccidioides immitis* IgG ELISA test using an oxidized and a non-oxidized antigen.

Figure 3 illustrates a *Aspergillus fumigatus* IgG ELISA test using an oxidized and a non-oxidized antigen.

5 Figure 4 illustrates a *Blastomyces dermatitidis* IgG ELISA test using an oxidized and a non-oxidized antigen.

Modes of Carrying Out the Invention

For clarity of disclosure, and not by way of limitation, the detailed description of
10 the invention is divided into the subsections that follow.

A. Definitions

Unless defined otherwise, all technical and scientific terms used herein have the same meaning as is commonly understood by one of ordinary skill in the art to which this invention belongs. All patents, applications, published applications and other publications referred to herein are incorporated by reference in their entirety. If a definition set forth in this section is contrary to or otherwise inconsistent with a definition set forth in the patents, applications, published applications and other publications that are herein incorporated by reference, the definition set forth in this section prevails over the definition that is incorporated herein by reference.

As used herein, “a” or “an” means “at least one” or “one or more.”

As used herein, “fungus” refers to a division of eucaryotic organisms that grow in irregular masses, without roots, stems, or leaves, and are devoid of chlorophyll or other pigments capable of photosynthesis. Each organism (thallus) is unicellular to filamentous, and possesses branched somatic structures (hyphae) surrounded by cell walls containing glucan or chitin or both, and containing true nuclei.

As used herein, “culture” refers to a growth of living cells or microorganisms, e.g., fungi, in a controlled artificial environment.

As used herein, “filtration” refers to a process of separating particulate matter from a fluid, such as air or liquid, by passing the fluid carrier through a medium that will not pass the particulates.

As used herein, “filtrate” refers to a fluid, *e.g.*, a liquid, that passes through the medium in filtration.

As used herein, “a fungal antigen as a culture filtrate” refers to a culture filtrate of a fungus containing the fungal antigenic substances. The culture filtrate can contain secreted extracellular fungal antigenic substances derived from supernatant of the fungal culture. Alternatively, culture filtrate can contain intracellular fungal antigenic substances derived from solution or suspension of disrupted fungal cells.

As used herein, “fungal antigen having (or has) not been purified” means that a particular fungal antigen (or a group of fungal antigens) has not been isolated or purified from other antigenic substances in the culture filtrate by any isolation or purification methods. However, the fungal antigen existing in the culture filtrate can be concentrated or diluted along with the concentration or dilution of the culture filtrate.

As used herein, “fungal antigen having not been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate” means that a particular fungal antigen (or a group of fungal antigens) has not been isolated or purified from other antigenic substances in the culture filtrate by ion exchange chromatography or isoelectric focusing. However, the fungal antigen existing in the culture filtrate can be concentrated or diluted along with the concentration or dilution of the culture filtrate, or can be isolated or purified from other antigenic substances in the culture filtrate by any isolation or purification methods other than ion exchange chromatography and isoelectric focusing.

As used herein, an “oxidizing agent (or oxidant)” refers to an agent that accepts electrons in an oxidation-reduction reaction.

As used herein the term “assessing (or assessed)” is intended to include quantitative and qualitative determination of the identity of a moiety, *e.g.*, an antigen, an antibody or an antigen-antibody complex, present in the sample, and also of obtaining an index, ratio, percentage, visual or other value indicative of the identity of a moiety in the sample. Assessment may be direct or indirect.

As used herein, “sample” refers to anything which may contain a fungal antigen or antibody to be tested by the present devices and/or methods. The sample may be a biological sample, such as a biological fluid or a biological tissue. Examples of biological fluids include urine, blood, plasma, serum, saliva, semen, stool, sputum,

cerebral spinal fluid, tears, mucus, amniotic fluid or the like. Biological tissues are aggregates of cells, usually of a particular kind together with their intercellular substance that form one of the structural materials of a human, animal, plant, bacterial, fungal or viral structure, including connective, epithelium, muscle and nerve tissues. Examples of 5 biological tissues also include organs, tumors, lymph nodes, arteries and individual cell(s). The sample may also be a fungal culture, *e.g.*, a fungal culture derived from a biological sample.

As used herein, a “liquid (fluid) sample” refers to a sample that naturally exists as a liquid or fluid, *e.g.*, a biological fluid. A “liquid sample” also refers to a sample that 10 naturally exists in a non-liquid status, *e.g.*, solid or gas, but is prepared as a liquid, fluid, solution or suspension containing the solid or gas sample material. For example, a liquid sample can encompass a liquid, fluid, solution or suspension containing a biological tissue.

15 **B. Methods for producing oxidized fungal antigen**

In one aspect, the present invention is directed to a method for producing a fungal antigen suitable for testing for an antibody to a fungus, the method comprising: a) 20 providing a fungal antigen as a culture filtrate, the fungal antigen having not been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate; and b) contacting the fungal antigen with an oxidizing agent to produce an oxidized fungal antigen suitable for testing for an antibody to the fungus.

The present methods can be used to produce an oxidized fungal antigen from any fungal genus or species. The oxidized antigens, while being reactive to antibodies to the fungus from which the fungal antigen is derived, may be reactive to antibodies to the 25 fungus of other related genus or species. Preferably, the methods will be controlled to produce oxidized fungal antigens that are able to distinguish among different fungal genus. More preferably, the methods will be controlled to produce oxidized fungal antigens that are able to distinguish among different fungal species within the same genus.

30 Fungi can be generally described as follows (Brock and Madigan, *Biology of Microorganisms (6th Ed.)*, Prentice Hall, Englewood Cliffs, New Jersey 07632 (1991),

pp. 817-818). In contrast to the algae, the fungi lack chlorophyll. Fungi can be differentiated from bacteria by the fact that fungal cells are usually much larger and contain a nucleus, vacuoles, and mitochondria, typical of eucaryotic cells. Although the fungi are a large and diverse group of eucaryotic microorganisms, three groups of fungi have major practical importance: the molds, yeasts, and mushrooms.

The habitats of fungi are quite diverse. Some are aquatic, living primarily in fresh water, and a few marine fungi are also known. Most fungi, however, have terrestrial habitats, in soil or on dead plant matter, and these types often play crucial roles in the mineralization of organic carbon in nature. A large number of fungi are parasites of terrestrial plants. Indeed, fungi cause the majority of economically significant diseases of crop plants (see Table 1). A few fungi are parasitic on animals, including humans, although in general fungi are less significant as animal pathogens than are bacteria and viruses.

15 **Table 1. Classification and major properties of fungi**

Group	Common name	Hyphae	Typical representatives	Type of sexual spore	Habitats	Common diseases
Ascomycetes	Sac fungi	Septate	Neurospora, Saccharomyces, Morchella (morels)	Ascospore	Soil, decaying plant material	Dutch elm, chestnut blight, ergot, rots
Basidiomycetes	Club fungi, mushrooms	Septate	Amanita (poisonous mushroom), Agaricus (edible mushroom)	Basidiospore	Soil, decaying plant material	Black stem, wheat rust, corn smut
Zygomycetes	Bread molds	Coenocytic	Mucor, Rhizopus (common bread mold)	Zygospore	Soil, decaying plant material	Food spoilage; rarely involved in parasitic disease
Oomycetes	Water molds	Coenocytic	Allomyces	Oospore	Aquatic	Potato blight, certain fish diseases

Deuteromyces	Fungi imperfecti	Septate	Penicillium, Aspergillus	None	Soil, decaying plant material, surfaces of animal bodies	Plant wilt, fungal infections of animals such as ringworm, athlete's foot, and other dermatomycoses
--------------	------------------	---------	--------------------------	------	--	---

All fungi are organotrophs. Lacking chlorophyll, they of course cannot photosynthesize, and the group also lacks lithotrophic forms. When compared to the bacteria, the fungi in general have fairly simple nutritional requirements, and their 5 metabolic and biosynthetic processes are not particularly diverse or unusual. It is in their morphological properties and in their sexual life cycles that the fungi exhibit considerable diversity; hence it is on the basis of these characteristics that the fungi are classified.

In a specific embodiment, the present method is used to produce an oxidized 10 fungal antigen from a mitosporic Trichocomaceae, an Onygenaceae or a mitosporic Onygenale. In another specific embodiment, the present method is used to produce an oxidized fungal antigen from *Aspergillus*, e.g., *Aspergillus fumigatus*, *Blastomyces (Ajellomyces)*, e.g., *Blastomyces dermatitidis*, *Coccidioides*, e.g., *Coccidioides immitis*, *Blastocystis*, e.g., *Blastocystis hominis*, *Histoplasma*, e.g., *Histoplasma capsulatum* and *Histoplasma duboifii*, *Candida*, e.g., *Candida albicans* and *Cryptococcus*, e.g., 15 *Cryptococcus neoformans*.

In still another specific embodiment, the present method is used to produce an oxidized fungal antigen from a fungal antigen not derived from *Aspergillus*, e.g., *Aspergillus fumigatus*, *Blastomyces (Ajellomyces)*, e.g., *Blastomyces dermatitidis*, *Coccidioides*, e.g., *Coccidioides immitis*, *Blastocystis*, e.g., *Blastocystis hominis*, *Histoplasma*, e.g., *Histoplasma capsulatum* and *Histoplasma duboifii*, *Candida*, e.g., *Candida albicans* and *Cryptococcus*, e.g., *Cryptococcus neoformans*. 20

The fungal antigen to be oxidized should not have been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate. However, the fungal 25 antigen to be oxidized can be concentrated or diluted along with the concentration or dilution of the culture filtrate, or can be isolated or purified from other antigenic substances in the culture filtrate by any isolation or purification methods other than ion

exchange chromatography and isoelectric focusing. For example, the fungal antigen can be isolated or purified by any chromatographic methods other than ion exchange chromatography or any electrophoresis other than isoelectric focusing, as well as other types of isolation or purification methods such as centrifugation or organic-aqueous phase separations. In a specific embodiment, the fungal antigen to be oxidized has not been purified from the culture filtrate prior to the oxidation step.

The present methods can comprise additional steps. For example, the present methods can further comprise a step of concentrating the culture filtrate. The concentrating step can be conducted prior to, concurrently with and/or subsequent to the oxidation step. Preferably, the concentrating step is conducted prior to the oxidation step. In another example, the present methods can further comprise a step of enzymatically deglycosylating the fungal antigen. Any suitable deglycosylating enzymes can be used, e.g., Endo H, PNGase F, *O*-glycanase, etc. The enzymatic deglycosylation step can be conducted prior to, concurrently with and/or subsequent to the oxidation step. Preferably, enzymatic deglycosylation step is conducted prior to the oxidation step.

Any suitable oxidizing agent can be used in the present methods. Exemplary oxidizing agents include hydrogen peroxide (H_2O_2), ozone (O_3), polyatomic oxygen O_7 , polyatomic oxygen O_8 , $NaIO_4$, potassium peroxymonosulfate (oxone) (Wozniak et al., *Bioorg. Med. Chem. Lett.*, 8(19):2641-6 (1998)), D,L-S-methylipoic acid methyl ester (Pan and Jordan, *Biochemistry*, 37(5):1357-64 (1998)), tertiary butyl hydroperoxide (Tarin et al., *Mol. Hum. Reprod.*, 2(12):895-901 (1996)), menadione (Santini et al., *Free Radic. Biol. Med.*, 20(7):915-24 (1996)), diamide (Bosin and Kasper, *J. Biochem. Toxicol.*, 7(3):139-45 (1992)), iodogen (Saha et al., *Int. J. Rad. Appl. Instrum.*, 16(4):431-3 (1989)), N-bromosuccinimide (Sinn et al., *Anal. Biochem.*, 170(1):186-92 (1988)), omeprazole (Im et al., *J. Biol. Chem.*, 260(8):4591-7 (1985)), and N-ethylmaleimide (Marzulli et al., *Boll. Soc. Ital. Biol. Sper.*, 61(1):121-7 (1985)).

In a specific embodiment, the oxidizing agent used in the present methods inactivates the cross-reactive carbohydrate moiety of the fungal antigen but does not disturb the structural integrity or antigenicity of the non-carbohydrate moieties. For example, when the fungal antigen to be oxidized is a glycosylated protein or peptide, the oxidizing agent inactivates the cross-reactive carbohydrate moiety of the fungal antigen

but does not disturb the structural integrity or antigenicity of the proteineous or peptidyl moiety.

In a preferred embodiment, the oxidizing agent used in the present methods is periodate (or periodic acid), which has a molecular formula H_5IO_6 . Suitable salts, esters or other derivatives of periodate that retain periodate's oxidizing activity, *e.g.*, sodium *meta*-periodate (Sigma). can also be used. Periodate or its salts, esters or other derivatives thereof can be produced by methods known in the art. For example, periodate can be prepared by electrolytic oxidation of iodic acid or from barium periodate and nitric acid (Willard, *Inorg. Syn.* 1, 172 (1939); Chemistry of periodic acid and periodates; H. Siebert, *Fortschr. Chem. Forsch.* 8, 470 (1967); Periodic acid and periodates in organic and bioorganic chemistry; A. J. Fatiadi, *Synthesis* 1974, 229; and G. Dryhurst, *Periodate Oxidation of Diol and Other Functional Groups* (Pergamon Press, New York, 1970)). Alternatively, Periodate or its salts, esters or other derivatives thereof can be obtained commercially, *e.g.*, from Sigma.

The oxidation step can be conducted under any suitable conditions depending on a number of factors such as the fungal antigen to be oxidized, the oxidizing agent to be used, the properties of the culture filtrate, the degree of concentration and/or dilution of the culture filtrate and the presence or absence of any enzymatic deglycosylation.

Generally, the concentration of the oxidizing agent can be in a wide range, *e.g.*, from 0.01 M to about 0.1 M if periodate is used. The oxidation step can be conducted at any suitable temperature, *e.g.*, from about 4°C to an ambient temperature, *e.g.*, from about 25°C to about 30°C. The incubation time can also have a wide range, *e.g.*, from about 1 hour to about 20 hours. When periodate is used, oxidation step is preferably conducted at about 25°C for about 1-3 hours.

In another specific embodiment, the present method can further comprise a step of attaching the oxidized fungal antigen to a surface of a device suitable for testing for an antibody to a fungus. Any suitable device can be used such as a microtiter plate, a glass slide, a nitrocellulose membrane, a latex bead, a cell, a test tube, a plastic bead, a colloidal gold particle, a colored particle, a magnetic bead and a quantum dot.

Preferably, a microtiter plate is used.

Fungal antigens, suitable for testing for an antibody to a fungus, which are produced by the present methods are also provided. Devices, suitable for testing for an antibody to a fungus, which are produced by the present methods are further provided. In a preferred embodiment, a device for simultaneously testing for a plurality of fungal antibodies is provided, the device comprising a plurality of fungal antigens produced by the present methods attached to a surface of the device suitable for testing for an antibody to a fungus, wherein the plurality of fungal antigens are attached to areas of the surface that are physically distinct from each other. Any suitable device can be used. In one example, the device comprises an immunoblot and wherein the plurality of fungal antigens are attached to separate stripes in the immunoblot. In another example, the device comprises an IFA well and wherein the plurality of fungal antigens are attached to different spots within the IFA well. In still another example, the device comprises an IFA slide and wherein the plurality of fungal antigens are attached to different wells of the IFA slide. Other exemplary devices include a microtiter plate, a glass slide, a nitrocellulose membrane, a latex bead, a cell, a test tube, a plastic bead, a colloidal gold particle, a colored particle, a magnetic bead and a quantum dot.

C. Methods of testing for an antibody to a fungus

In another aspect, the present invention is directed to a method of testing for an antibody to a fungus in a sample, the method comprising: a) producing a fungal antigen suitable for testing for an antibody to a fungus, comprising providing a fungal antigen as a culture filtrate, the fungal antigen having not been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate, and contacting the fungal antigen with an oxidizing agent to produce an oxidized fungal antigen suitable for testing for an antibody to the fungus; b) contacting a sample suspected of containing an antibody to a fungus with the oxidized fungal antigen produced in step a) under suitable conditions to allow binding of the antibody, if present in the sample, to the oxidized fungal antigen; and c) assessing binding between the antibody and the oxidized fungal antigen to determine the presence and/or amount of the antibody in the sample.

The present methods can be used to test for an antibody to any fungus in a sample, including any fungus described in the above Section B. Preferably, the antibody

to be tested is to a pathogenic fungus such as *Aspergillus*, e.g., *Aspergillus fumigatus*, *Blastomyces* (*Ajellomyces*), e.g., *Blastomyces dermatitidis*, *Coccidioides*, e.g., *Coccidioides immitis*, *Blastocystis*, e.g., *Blastocystis hominis*, *Histoplasma*, e.g., *Histoplasma capsulatum* and *Candida*, e.g., *Candida albicans* and *Candida neoformans*.

5 The present methods can be used to test any sample. Preferably, the sample to be tested is a clinical sample including human and veterinary clinical samples. More preferably, the sample to be tested is a human clinical sample.

10 The oxidized fungal antigens to be used in the present test methods are produced as described in the above Section B. Preferably, the fungal antigen to be oxidized has not been purified from the culture filtrate. Also, preferably, the oxidizing agent to be used is periodate.

15 The present test methods can be conducted in liquid or solid phase. In specific embodiment, the present test method is conducted in a solid phase the method further comprises a step of attaching the oxidized fungal antigen to a surface of a device suitable for testing for an antibody to a fungus before contacting the antigen with the sample.

20 The binding between the antibody and the oxidized fungal antigen can be assessed by any suitable assay formats. For example, the binding between the antibody and the oxidized fungal antigen can be assessed by a sandwich or competitive assay format. In another example, the binding between the antibody and the oxidized fungal antigen can be assessed by an enzyme-linked immunosorbent assay (ELISA), immunoblotting, immunoprecipitation, radioimmunoassay (RIA), immunostaining, latex agglutination, indirect hemagglutination assay (IHA), complement fixation, indirect immunofluorescent assay (IFA), nephelometry, flow cytometry assay, chemiluminescence assay, lateral flow immunoassay, u-capture assay, inhibition assay or avidity assay (Manual of Clinical 25 Laboratory Immunology, NR Rose, E Conway de Macario, JD Folds, HC Lane, RM Nakamura, eds. ASM Press, Washington DC, 1997). Preferably, the binding between the antibody and the oxidized fungal antigen is assessed by an ELISA format.

30 The present assay method can be conducted for prognosis, diagnosis and/or monitoring treatment of pathogenic fungal infection, e.g., histoplasmosis caused by *Histoplasma capsulatum* infection, blastomycosis caused by *Blastomyces dermatitidis* infection, candidiasis caused by *Candida albicans* or *Candida neoformans* infection, etc.

D. Methods for producing an antibody to an fungal antigen

In still another aspect, the present invention is directed to a method for producing an antibody to a fungal antigen, the method comprising: a) producing a fungal antigen comprising providing a fungal antigen as a culture filtrate, the fungal antigen having not been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate, and contacting the fungal antigen with an oxidizing agent to produce an oxidized fungal antigen; b) delivering, to a vertebrate or tissue culture, the oxidized fungal antigen, in an amount sufficient to induce detectable production of an antibody to the antigen; and c) recovering the antibody from the vertebrate or tissue culture.

The oxidized fungal antigens to be used in the present antibody producing methods are produced as described in the above Section B. Preferably, the fungal antigen to be oxidized has not been purified from the culture filtrate. Also, preferably, the oxidizing agent to be used is periodate.

Any suitable non-human mammal can be used in the present methods. For example, mouse, rabbit and goat can be used.

The antibody to a fungal antigen, *e.g.*, a polyclonal antiserum, which is produced by the present method is also provided.

In yet another aspect, the present invention is directed to a method for producing a monoclonal antibody to an fungal antigen, the method comprising: a) producing a fungal antigen comprising providing a fungal antigen as a culture filtrate, the fungal antigen having not been purified by ion exchange chromatography or isoelectric focusing from the culture filtrate, and contacting the fungal antigen with an oxidizing agent to produce an oxidized fungal antigen; b) delivering, to a vertebrate or tissue culture, the oxidized fungal antigen, in an amount sufficient to induce detectable production of an antibody to the antigen; c) removing at least a portion of antibody-producing cells from the vertebrate or tissue culture; d) immortalizing the removed antibody-producing cells; e) propagating the immortalized antibody-producing cells; and f) harvesting monoclonal antibody produced by the immortalized antibody-producing cells. Methods of testing for a fungal antigen in a sample using the produced monoclonal antibodies are also provided.

The oxidized fungal antigens to be used in the present antibody producing methods are produced as described in the above Section B. Preferably, the fungal antigen to be oxidized has not been purified from the culture filtrate. Also, preferably, the oxidizing agent to be used is periodate.

5 Any suitable non-human mammal can be used in the present methods. For example, mouse, rabbit and goat can be used.

A monoclonal antibody to a fungal antigen, which is produced by the present method is also provided.

10 A hybridoma capable of producing a monoclonal antibody to a fungal antigen, which is produced by steps a)-d) of the present method is further provided.

The oxidized fungal antigens can be delivered to a vertebrate or non-human mammal by any methods known in the art (*See e.g.*, Coligan et al. (Ed.), Current Protocols in Immunology, 2.II Production of Antibodies, John Wiley & Sons, Inc. (2000)).

15 The oxidized fungal antigens can be delivered to the interstitial space of tissues of the animal body, including those of muscle, skin, brain, lung, liver, spleen, bone marrow, thymus, heart lymph, blood, bone, cartilage, pancreas, kidney, gall bladder, stomach, intestine, testis, ovary, uterus, rectum, nervous system, eye, gland, and connective tissue. Interstitial space of the tissues comprises the intercellular, fluid, mucopolysaccharide matrix among the reticular fibers or organ tissues, elastic fibers in the walls of vessels or chambers, collagen fibers of fibrous tissues, or that same matrix within connective tissue ensheathing muscle cells or in the lacunae of bone. It is similarly the space occupied by the plasma of the circulation of the lymph fluid of the lymphatic channels.

20 The oxidized fungal antigens can be conveniently delivered by injection into the tissues comprising these cells. They are preferably delivered to persistent, non-dividing cells which are differentiated, although delivery can be achieved in non-differentiated or less completely differentiated cells, such as, for example, stem cells of blood or skin fibroblasts.

25 In a specific embodiment, the oxidized fungal antigen is delivered directly to a tissue of the animals. Preferably, the oxidized fungal antigen is delivered directly to muscle, skin or mucous membrane. In one example, the oxidized fungal antigen can be

delivered directly to a tissue of the animal by injection, by gene gun technology or by lipid mediated delivery technology. The injection can be conducted via a needle or other injection devices. The gene gun technology is disclosed in U.S. Patent No. 5,302,509 and the lipid mediated delivery technology is disclosed in U.S. Patent No. 5,703,055, the

5 contents of which are incorporated herein by reference.

In still another specific embodiment, the oxidized fungal antigen is delivered to a cell of the animal and the cell containing the oxidized fungal antigen is delivered to a suitable tissue of the animal. Preferably, the oxidized fungal antigen is delivered to a blood cell of an animal. More preferably, the oxidized fungal antigen is delivered to a

10 spleen B cell of an animal.

The anti-fungal antibodies, whether polyclonal or monoclonal, that are produced by the above methods can be used in testing for a fungal antigen in a sample. In one specific embodiment, the present invention is directed to a method of testing for a fungal antigen in a sample, the method comprising: a) providing an anti-fungal polyclonal

15 antiserum produced by the above antibody producing method; b) contacting a sample suspected of containing a fungal antigen with the polyclonal antiserum under suitable conditions to allow binding of the fungal antigen, if present in the sample, to the polyclonal antiserum; and c) assessing binding between the fungal antigen and the polyclonal antiserum to determine the presence and/or amount of the fungal antigen in the sample.

In another specific embodiment, the present invention is directed to a method of testing for a fungal antigen in a sample, the method comprising: a) providing an anti-fungal antigen monoclonal antibody produced by the above antibody producing method;

20 b) contacting a sample suspected of containing a fungal antigen with the monoclonal antibody under suitable conditions to allow binding of the fungal antigen, if present in the sample, to the monoclonal antibody; and c) assessing binding between the fungal antigen and the monoclonal antibody to determine the presence and/or amount of the fungal antigen in the sample.

30 **E. Exemplary embodiments**

In exemplary embodiments of the present invention, oxidized preparations of crude fungal antigens are used in an ELISA format, resulting in the sensitive and specific detection of fungal antibodies. The fungal antigens are oxidized during a short (2-hour) incubation with periodate, then further diluted with coating buffer and allowed to adsorb to polystyrene microtiter wells. ELISA performance using these coated microtiter wells is complete within 4 hours, thus providing sensitive and specific results more quickly and easily than ID and western blot assays.

Preparation of antigen-coated microtiter wells

A fungal antigen preparation, designed for use in ID assays, is purchased from a commercial source. The protein concentration of the antigen preparation is determined by protein-dye binding (Bradford) assay. The fungal antigen is then diluted to a final protein concentration of about 10-50 micrograms per milliliter (ug/mL) in 0.05M sodium acetate buffer (adjusted to pH 5.5 using 10% acetic acid) containing about 0.02M-0.10M sodium meta-periodate. After a short incubation period (1-3 hours) at room temperature, during which hydroxyl groups of the fungal antigen carbohydrate moieties are cleaved and converted to aldehydes (i.e., become oxidized), the oxidized fungal antigen preparation is diluted approximately 20-fold with phosphate buffered saline (PBS), pH 7.4. This diluted fungal antigen is then added to polystyrene microtiter wells at a volume of 0.1 mL per well. The microtiter wells are covered with plastic adhesive tape and incubated in the refrigerator (4-8°C) for 14-18 hours. During this incubation, fungal antigen protein moieties attach to the polystyrene through hydrophobic interactions. After incubation, the solution is discarded from the microtiter wells, now bearing attached oxidized fungal antigen. All microtiter wells then receive 0.2 mL of PBS containing a blocking agent, such as bovine serum albumin (BSA), and a stabilizing agent, such as trehalose. The blocking agent attaches to any polystyrene binding sites that are left exposed after the incubation with oxidized fungal antigen. This step reduces the nonspecific attachment of antibody molecules not directed toward the specific antigens, and also reduces nonspecific adsorption of the conjugate during the color generation step. After 2 hours at room temperature, the blocking solution is discarded from the microtiter wells. The microtiter plates are air-dried at room temperature, then

placed inside a plastic pouch, and the pouch heat-sealed. Sealed pouches containing microtiter plates coated with fungal antigens are stored in the refrigerator until use in the ELISA.

5

ELISA performance

Antigen-coated microtiter well plates are removed from the refrigerator and allowed to warm to room temperature. The wells are then filled with PBS containing a mild detergent such as Tween 20 (PBST), and after 5 minutes this PBST is discarded. This wash step ensures that all antigen and blocking solution proteins have been removed; in addition, the small amount of detergent helps to reduce bubble formation in subsequent steps. Serum samples are diluted in PBST that contains the same protein that is used as a blocking agent (such as BSA). About 0.1 mL of each diluted serum is added to an assigned microtiter well. The microtiter plate is then covered with adhesive tape and incubated about 1.5 hours at room temperature. During this incubation period, any serum antibodies directed to the antigen that is attached to the microtiter well will bind to the antigen. At the end of this incubation, the liquid material is discarded from all the wells. About 0.25 mL of PBST is then added to each well, and this fluid is then discarded; this wash step is repeated three times. This washing procedure ensures that non-attached serum components are removed from the microtiter well. About 0.1 mL of a conjugate diluted in PBST is then added to each well. This conjugate is a commercially-available, species-specific, anti-immunoglobulin with an enzyme chemically linked (conjugated) to it. The enzyme used is selected based on its ability to react with a substrate to produce a colored compound. A typical conjugate would be goat anti-human IgG tagged with horseradish peroxidase. This conjugate will bind to any IgG antibodies in the serum that attached to the antigen attached to the microtiter well. Following another series of washes to remove any unbound conjugate, enzyme substrate is added to each microtiter well. A typical substrate would be tetra-methylbenzidine (TMB) and hydrogen peroxide; in the presence of horseradish peroxidase, this substrate produces a bluish purple color. Thus, color only develops when a "sandwich" consisting of plastic-bound antigen, serum-derived antibody, and enzyme-tagged conjugate has been built up; if the serum contains no antigen-specific antibody molecules, then the enzyme-

tagged conjugate will not bind, and no color will develop because the enzyme is absent from the system. After a defined amount of time, a solution of dilute hydrochloric acid, or some other acid, is added to lower the pH of the solution in the wells. At a lower pH, the bluish purple color changes to a yellow color. The absorbance (also called optical density) of the yellow color is then quantitated using a spectrophotometric instrument. The microtiter well plate is placed inside the instrument, and the optical density for each well (i.e., each sample) is measured at a light wavelength of 450 nanometers. The optical density value is directly proportional to the intensity of the yellow color, which in turn is directly proportional to the amount of antigen-specific IgG antibody bound to the well.

10

EXAMPLE 1: Detection of IgG antibodies recognizing *Histoplasma capsulatum*

Histoplasma capsulatum antigen (culture filtrate containing both H and M antigens) intended for use in ID assays was purchased from a commercial source, and adjusted to 20 μ g/mL in sodium acetate buffer containing 0.05M periodate (hereafter referred to as oxidized antigen). For comparative purposes, a duplicate antigen suspension in sodium acetate buffer without periodate was prepared (hereafter referred to as nonoxidized antigen). After 2 hours at room temperature, the antigen preparations were diluted 20-fold in PBS to a final concentration of 1 μ g/mL, then added to microtiter wells (0.1 mL per well). Thus, some microtiter wells received oxidized Histoplasma antigen, whereas other wells received nonoxidized Histoplasma antigen. After an overnight incubation in the refrigerator, the Histoplasma antigen preparations (oxidized or nonoxidized) were discarded from the wells, and PBS containing 0.1% (w/v) BSA and 5% (w/v) trehalose was added to the wells as blocking agent. After 2 hours at room temperature, the blocking agent was discarded, and the plates were then air-dried at room temperature for 2 hours. After pouching, the plates were stored in the refrigerator.

25
30
35
40
45
50
55
60
65
70
75
80
85
90
95
100
105
110
115
120
125
130
135
140
145
150
155
160
165
170
175
180
185
190
195
200
205
210
215
220
225
230
235
240
245
250
255
260
265
270
275
280
285
290
295
300
305
310
315
320
325
330
335
340
345
350
355
360
365
370
375
380
385
390
395
400
405
410
415
420
425
430
435
440
445
450
455
460
465
470
475
480
485
490
495
500
505
510
515
520
525
530
535
540
545
550
555
560
565
570
575
580
585
590
595
600
605
610
615
620
625
630
635
640
645
650
655
660
665
670
675
680
685
690
695
700
705
710
715
720
725
730
735
740
745
750
755
760
765
770
775
780
785
790
795
800
805
810
815
820
825
830
835
840
845
850
855
860
865
870
875
880
885
890
895
900
905
910
915
920
925
930
935
940
945
950
955
960
965
970
975
980
985
990
995
1000
1005
1010
1015
1020
1025
1030
1035
1040
1045
1050
1055
1060
1065
1070
1075
1080
1085
1090
1095
1100
1105
1110
1115
1120
1125
1130
1135
1140
1145
1150
1155
1160
1165
1170
1175
1180
1185
1190
1195
1200
1205
1210
1215
1220
1225
1230
1235
1240
1245
1250
1255
1260
1265
1270
1275
1280
1285
1290
1295
1300
1305
1310
1315
1320
1325
1330
1335
1340
1345
1350
1355
1360
1365
1370
1375
1380
1385
1390
1395
1400
1405
1410
1415
1420
1425
1430
1435
1440
1445
1450
1455
1460
1465
1470
1475
1480
1485
1490
1495
1500
1505
1510
1515
1520
1525
1530
1535
1540
1545
1550
1555
1560
1565
1570
1575
1580
1585
1590
1595
1600
1605
1610
1615
1620
1625
1630
1635
1640
1645
1650
1655
1660
1665
1670
1675
1680
1685
1690
1695
1700
1705
1710
1715
1720
1725
1730
1735
1740
1745
1750
1755
1760
1765
1770
1775
1780
1785
1790
1795
1800
1805
1810
1815
1820
1825
1830
1835
1840
1845
1850
1855
1860
1865
1870
1875
1880
1885
1890
1895
1900
1905
1910
1915
1920
1925
1930
1935
1940
1945
1950
1955
1960
1965
1970
1975
1980
1985
1990
1995
2000
2005
2010
2015
2020
2025
2030
2035
2040
2045
2050
2055
2060
2065
2070
2075
2080
2085
2090
2095
2100
2105
2110
2115
2120
2125
2130
2135
2140
2145
2150
2155
2160
2165
2170
2175
2180
2185
2190
2195
2200
2205
2210
2215
2220
2225
2230
2235
2240
2245
2250
2255
2260
2265
2270
2275
2280
2285
2290
2295
2300
2305
2310
2315
2320
2325
2330
2335
2340
2345
2350
2355
2360
2365
2370
2375
2380
2385
2390
2395
2400
2405
2410
2415
2420
2425
2430
2435
2440
2445
2450
2455
2460
2465
2470
2475
2480
2485
2490
2495
2500
2505
2510
2515
2520
2525
2530
2535
2540
2545
2550
2555
2560
2565
2570
2575
2580
2585
2590
2595
2600
2605
2610
2615
2620
2625
2630
2635
2640
2645
2650
2655
2660
2665
2670
2675
2680
2685
2690
2695
2700
2705
2710
2715
2720
2725
2730
2735
2740
2745
2750
2755
2760
2765
2770
2775
2780
2785
2790
2795
2800
2805
2810
2815
2820
2825
2830
2835
2840
2845
2850
2855
2860
2865
2870
2875
2880
2885
2890
2895
2900
2905
2910
2915
2920
2925
2930
2935
2940
2945
2950
2955
2960
2965
2970
2975
2980
2985
2990
2995
3000
3005
3010
3015
3020
3025
3030
3035
3040
3045
3050
3055
3060
3065
3070
3075
3080
3085
3090
3095
3100
3105
3110
3115
3120
3125
3130
3135
3140
3145
3150
3155
3160
3165
3170
3175
3180
3185
3190
3195
3200
3205
3210
3215
3220
3225
3230
3235
3240
3245
3250
3255
3260
3265
3270
3275
3280
3285
3290
3295
3300
3305
3310
3315
3320
3325
3330
3335
3340
3345
3350
3355
3360
3365
3370
3375
3380
3385
3390
3395
3400
3405
3410
3415
3420
3425
3430
3435
3440
3445
3450
3455
3460
3465
3470
3475
3480
3485
3490
3495
3500
3505
3510
3515
3520
3525
3530
3535
3540
3545
3550
3555
3560
3565
3570
3575
3580
3585
3590
3595
3600
3605
3610
3615
3620
3625
3630
3635
3640
3645
3650
3655
3660
3665
3670
3675
3680
3685
3690
3695
3700
3705
3710
3715
3720
3725
3730
3735
3740
3745
3750
3755
3760
3765
3770
3775
3780
3785
3790
3795
3800
3805
3810
3815
3820
3825
3830
3835
3840
3845
3850
3855
3860
3865
3870
3875
3880
3885
3890
3895
3900
3905
3910
3915
3920
3925
3930
3935
3940
3945
3950
3955
3960
3965
3970
3975
3980
3985
3990
3995
4000
4005
4010
4015
4020
4025
4030
4035
4040
4045
4050
4055
4060
4065
4070
4075
4080
4085
4090
4095
4100
4105
4110
4115
4120
4125
4130
4135
4140
4145
4150
4155
4160
4165
4170
4175
4180
4185
4190
4195
4200
4205
4210
4215
4220
4225
4230
4235
4240
4245
4250
4255
4260
4265
4270
4275
4280
4285
4290
4295
4300
4305
4310
4315
4320
4325
4330
4335
4340
4345
4350
4355
4360
4365
4370
4375
4380
4385
4390
4395
4400
4405
4410
4415
4420
4425
4430
4435
4440
4445
4450
4455
4460
4465
4470
4475
4480
4485
4490
4495
4500
4505
4510
4515
4520
4525
4530
4535
4540
4545
4550
4555
4560
4565
4570
4575
4580
4585
4590
4595
4600
4605
4610
4615
4620
4625
4630
4635
4640
4645
4650
4655
4660
4665
4670
4675
4680
4685
4690
4695
4700
4705
4710
4715
4720
4725
4730
4735
4740
4745
4750
4755
4760
4765
4770
4775
4780
4785
4790
4795
4800
4805
4810
4815
4820
4825
4830
4835
4840
4845
4850
4855
4860
4865
4870
4875
4880
4885
4890
4895
4900
4905
4910
4915
4920
4925
4930
4935
4940
4945
4950
4955
4960
4965
4970
4975
4980
4985
4990
4995
5000
5005
5010
5015
5020
5025
5030
5035
5040
5045
5050
5055
5060
5065
5070
5075
5080
5085
5090
5095
5100
5105
5110
5115
5120
5125
5130
5135
5140
5145
5150
5155
5160
5165
5170
5175
5180
5185
5190
5195
5200
5205
5210
5215
5220
5225
5230
5235
5240
5245
5250
5255
5260
5265
5270
5275
5280
5285
5290
5295
5300
5305
5310
5315
5320
5325
5330
5335
5340
5345
5350
5355
5360
5365
5370
5375
5380
5385
5390
5395
5400
5405
5410
5415
5420
5425
5430
5435
5440
5445
5450
5455
5460
5465
5470
5475
5480
5485
5490
5495
5500
5505
5510
5515
5520
5525
5530
5535
5540
5545
5550
5555
5560
5565
5570
5575
5580
5585
5590
5595
5600
5605
5610
5615
5620
5625
5630
5635
5640
5645
5650
5655
5660
5665
5670
5675
5680
5685
5690
5695
5700
5705
5710
5715
5720
5725
5730
5735
5740
5745
5750
5755
5760
5765
5770
5775
5780
5785
5790
5795
5800
5805
5810
5815
5820
5825
5830
5835
5840
5845
5850
5855
5860
5865
5870
5875
5880
5885
5890
5895
5900
5905
5910
5915
5920
5925
5930
5935
5940
5945
5950
5955
5960
5965
5970
5975
5980
5985
5990
5995
6000
6005
6010
6015
6020
6025
6030
6035
6040
6045
6050
6055
6060
6065
6070
6075
6080
6085
6090
6095
6100
6105
6110
6115
6120
6125
6130
6135
6140
6145
6150
6155
6160
6165
6170
6175
6180
6185
6190
6195
6200
6205
6210
6215
6220
6225
6230
6235
6240
6245
6250
6255
6260
6265
6270
6275
6280
6285
6290
6295
6300
6305
6310
6315
6320
6325
6330
6335
6340
6345
6350
6355
6360
6365
6370
6375
6380
6385
6390
6395
6400
6405
6410
6415
6420
6425
6430
6435
6440
6445
6450
6455
6460
6465
6470
6475
6480
6485
6490
6495
6500
6505
6510
6515
6520
6525
6530
6535
6540
6545
6550
6555
6560
6565
6570
6575
6580
6585
6590
6595
6600
6605
6610
6615
6620
6625
6630
6635
6640
6645
6650
6655
6660
6665
6670
6675
6680
6685
6690
6695
6700
6705
6710
6715
6720
6725
6730
6735
6740
6745
6750
6755
6760
6765
6770
6775
6780
6785
6790
6795
6800
6805
6810
6815
6820
6825
6830
6835
6840
6845
6850
6855
6860
6865
6870
6875
6880
6885
6890
6895
6900
6905
6910
6915
6920
6925
6930
6935
6940
6945
6950
6955
6960
6965
6970
6975
6980
6985
6990
6995
7000
7005
7010
7015
7020
7025
7030
7035
7040
7045
7050
7055
7060
7065
7070
7075
7080
7085
7090
7095
7100
7105
7110
7115
7120
7125
7130
7135
7140
7145
7150
7155
7160
7165
7170
7175
7180
7185
7190
7195
7200
7205
7210
7215
7220
7225
7230
7235
7240
7245
7250
7255
7260
7265
7270
7275
7280
7285
7290
7295
7300
7305
7310
7315
7320
7325
7330
7335
7340
7345
7350
7355
7360
7365
7370
7375
7380
7385
7390
7395
7400
7405
7410
7415
7420
7425
7430
7435
7440
7445
7450
7455
7460
7465
7470
7475
7480
7485
7490
7495
7500
7505
7510
7515
7520
7525
7530
7535
7540
7545
7550
7555
7560
7565
7570
7575
7580
7585
7590
7595
7600
7605
7610
7615
7620
7625
7630
7635
7640
7645
7650
7655
7660
7665
7670
7675
7680
7685
7690
7695
7700
7705
7710
7715
7720
7725
7730
7735
7740
7745
7750
7755
7760
7765
7770
7775
7780
7785
7790
7795
7800
7805
7810
7815
7820
7825
7830
7835
7840
7845
7850
7855
7860
7865
7870
7875
7880
7885
7890
7895
7900
7905
7910
7915
7920
7925
7930
7935
7940
7945
7950
7955
7960
7965
7970
7975
7980
7985
7990
7995
8000
8005
8010
8015
8020
8025
8030
8035
8040
8045
8050
8055
8060
8065
8070
8075
8080
8085
8090
8095
8100
8105
8110
8115
8120
8125
8130
8135
8140
8145
8150
8155
8160
8165
8170
8175
8180
8185
8190
8195
8200
8205
8210
8215
8220
8225
8230
8235
8240
8245
8250
8255
8260
8265
8270
8275
8280
8285
8290
8295
8300
8305
8310
8315
8320
8325
8330
8335
8340
8345
8350
8355
8360
8365
8370
8375
8380
8385
8390
8395
8400
8405
8410
8415
8420
8425
8430
8435
8440
8445
8450
8455
8460
8465
8470
8475
8480
8485
8490
8495
8500
8505
8510
8515
8520
8525
8530
8535
8540
8545
8550
8555
8560
8565
8570
8575
8580
8585
8590
8595
8600
8605
8610
8615
8620
8625
8630
8635
8640
8645
8650
8655
8660
8665
8670
8675
8680
8685
8690
8695
8700
8705
8710
8715
8720
8725
8730
8735
8740
8745
8750
8755
8760
8765
8770
8775
8780
8785
8790
8795
8800
8805
8810
8815
8820
8825
8830
8835
8840
8845
8850
8855
8860
8865
8870
8875
8880
8885
8890
8895
8900
8905
8910
8915
8920
8925
8930
8935
8940
8945
8950
8955
8960
8965
8970
8975
8980
8985
8990
8995
9000
9005
9010
9015
9020
9025
9030
9035
9040
9045
9050
9055
9060
9065
9070
9075
9080
9085
9090
9095
9100
9105
9110
9115
9120
9125
9130
9135
9140
9145

serum samples were then added to duplicate microtiter wells (one containing oxidized Histoplasma antigen, the other containing nonoxidized Histoplasma antigen). The microtiter wells were covered with adhesive tape and incubated at room temperature for 1.5 hours. Following 4 washes with PBST, each well then received 0.1 mL of

5 horseradish peroxidase-conjugated goat anti-human IgG (diluted 1:10000 in PBST), and incubation was continued for 1 hour at room temperature. Following another series of 4 washes, each well then received 0.1 mL of tetra-methylbenzidine; after 10 minutes at room temperature, the reaction was stopped by the addition of 0.1N hydrochloric acid. The optical density was measured at 450 nM using a spectrophotometric ELISA reader.

10 Figure 1 presents the optical density (OD) values obtained following incubation of ID-positive and ID-negative serum samples in microtiter wells coated with nonoxidized and oxidized Histoplasma antigen. When using nonoxidized Histoplasma antigen as substrate, the distribution of OD values for the ID-negative serum group was similar to that observed for the ID-positive serum group. In contrast, when using the

15 oxidized Histoplasma antigen as substrate, the distribution of values for the ID-negative serum group was markedly different from the distribution of values for the ID-positive serum group. When using the oxidized antigen, 16 of 16 (100%) ID-positive serum samples exhibited OD values >1.45, whereas only 4 of 64 (6%) ID-negative serum samples exhibited OD values >1.45. Expressed another way, the Histoplasma IgG assay using oxidized Histoplasma antigen as substrate was 100% sensitive and 94% specific

20 (see Table 2) using an OD cutoff value of 1.45.

25 Table 2. Specificities of fungal IgG ELISAs using oxidized versus nonoxidized antigens under conditions giving the same sensitivities

ELISA	Antigen type	Cutoff OD	Sensitivity (%)	Specificity (%)
Histoplasma IgG	Oxidized	1.45	100	94
	Nonoxidized	0.50	100	25
Coccidioides IgG	Oxidized	1.00	90	97
	Nonoxidized	0.22	90	28
Aspergillus IgG	Oxidized	.85	89	95
	Nonoxidized	0.23	89	65
Blastomyces IgG	Oxidized	0.40	100	100
	Nonoxidized	0.50	100	87

For comparative purposes, we determined the specificity of the Histoplasma IgG assay using nonoxidized antigen under conditions where the sensitivity using nonoxidized antigen matched that seen using oxidized antigen (i.e., 100%). Thus, the cutoff OD giving 100% sensitivity when using nonoxidized antigen was 0.50 (Figure 1).

5 Based on this cutoff OD value, the specificity of the Histoplasma IgG assay using nonoxidized antigen was only 25% (see Table 2). Thus, under conditions providing 100% sensitivity, oxidation of Histoplasma antigen profoundly improved the specificity of an ELISA for detecting Histoplasma IgG.

10 **EXAMPLE 2: Detection of IgG antibodies recognizing *Coccidioides immitis***

15 *Coccidioides immitis* antigen (culture filtrate containing F antigen) intended for use in ID assays was purchased from a commercial source. Microtiter well preparation and assay performance was the same as described in Example 1 for Histoplasma, with the exception that conjugate was used at a final dilution of 1:5000 rather than 1:10000. The serum panel evaluated for Coccidioides IgG consisted of 80 samples previously tested for Coccidioides antibodies by ID (41 positive, 39 negative).

20 Figure 2 presents the optical density (OD) values obtained following incubation of ID-positive and ID-negative serum samples in microtiter wells coated with nonoxidized and oxidized Coccidioides antigen. When using nonoxidized Coccidioides antigen as substrate, the distribution of OD values for the ID-negative serum group was similar to that observed for the ID-positive serum group. In contrast, when using the oxidized Coccidioides antigen as substrate, the distribution of values for the ID-negative serum group was markedly different from the distribution of values for the ID-positive serum group. When using the oxidized antigen, 37 of 41 (90%) ID-positive serum samples exhibited OD values >1.00, whereas only 1 of 39 (3%) ID-negative serum samples exhibited OD values >1.00. Expressed another way, the Coccidioides IgG assay using oxidized Coccidioides antigen as substrate was 90% sensitive and 97% specific (see Table 2) using an OD cutoff value of 1.00.

25 For comparative purposes, we determined the specificity of the Coccidioides IgG assay using nonoxidized antigen under conditions where the sensitivity using nonoxidized antigen matched that seen using oxidized antigen (i.e., 90%). Thus, the

cutoff OD giving 90% sensitivity when using nonoxidized antigen was 0.22 (Figure 2). Based on this cutoff OD value, the specificity of the *Coccidioides* IgG assay using nonoxidized antigen was only 28% (see Table 2). Thus, under conditions providing 90% sensitivity, oxidation of *Coccidioides* antigen profoundly improved the specificity of an 5 ELISA for detecting *Coccidioides* IgG.

EXAMPLE 3: Detection of IgG antibodies recognizing *Aspergillus fumigatus*

Aspergillus fumigatus antigen (culture filtrate) intended for use in ID assays was purchased from a commercial source. Microtiter well preparation and assay performance 10 was the same as described in Example 1 for *Histoplasma*, with the exception that conjugate was used at a final dilution of 1:5000 rather than 1:10000. The serum panel evaluated for *Aspergillus* IgG consisted of 29 samples previously tested for *Aspergillus* antibodies by ID (9 positive, 20 negative).

Figure 3 presents the optical density (OD) values obtained following incubation 15 of ID-positive and ID-negative serum samples in microtiter wells coated with nonoxidized and oxidized *Aspergillus* antigen. For both nonoxidized and oxidized antigen, the overall distribution of OD values for the ID-negative group was visually distinct from the distribution of OD values for the ID-positive group. However, the range of difference between the ID-negative distribution and the ID-positive distribution was 20 more marked when using oxidized *Aspergillus* antigen. When using the oxidized antigen, 8 of 9 (89%) ID-positive serum samples exhibited OD values >0.85 , whereas only 1 of 20 (5%) ID-negative serum samples exhibited OD values >0.85 . Expressed another way, the *Aspergillus* IgG assay using oxidized *Aspergillus* antigen as substrate was 89% sensitive and 95% specific (see Table 2) using an OD cutoff value of 0.85.

25 For comparative purposes, we determined the specificity of the *Aspergillus* IgG assay using nonoxidized antigen under conditions where the sensitivity using nonoxidized antigen matched that seen using oxidized antigen (i.e., 89%). Thus, the cutoff OD giving 89% sensitivity when using nonoxidized antigen was 0.23 (Figure 3). Based on this cutoff OD value, the specificity of the *Aspergillus* IgG assay using 30 nonoxidized antigen was 65% (see Table 2). Thus, under conditions providing 89%

sensitivity, oxidation of Aspergillus antigen markedly improved the specificity of an ELISA for detecting Aspergillus IgG.

EXAMPLE 4: Detection of IgG antibodies recognizing *Blastomyces dermatitidis*

5 *Blastomyces dermatitidis* antigen (culture filtrate containing A antigen) intended for use in ID assays was purchased from a commercial source. Microtiter well preparation and assay performance was exactly the same as described in Example 1 for Histoplasma. The serum panel evaluated for Blastomyces IgG consisted of 16 samples previously tested for Blastomyces antibodies by ID (1 positive, 15 negative).

10 Figure 4 presents the optical density (OD) values obtained following incubation of ID-positive and ID-negative serum samples in microtiter wells coated with nonoxidized and oxidized Blastomyces antigen. When using the oxidized antigen, the single ID-positive serum sample exhibited an OD value >0.4 , whereas only 0 of 15 (0%) ID-negative serum samples exhibited OD values >0.40 . Expressed another way, the 15 Blastomyces IgG assay using oxidized Basltomyces antigen as substrate was 100% sensitive and 100% specific (see Table 2) using an OD cutoff value of 0.40.

For comparative purposes, we determined the specificity of the Blastomyces IgG assay using nonoxidized antigen under conditions where the sensitivity using nonoxidized antigen matched that seen using oxidized antigen (i.e., 100%). Thus, the 20 cutoff OD giving 100% sensitivity when using nonoxidized antigen was 0.50 (Figure 4). Based on this cutoff OD value, the specificity of the Blastomyces IgG assay using nonoxidized antigen was 87% (see Table 2). Thus, under conditions providing 100% sensitivity, oxidation of Blastomyces antigen improved the specificity of an ELISA for detecting Blastomyces IgG.

25 The above examples are included for illustrative purposes only and are not intended to limit the scope of the invention. Many variations to those described above are possible. Since modifications and variations to the examples described above will be apparent to those of skill in this art, it is intended that this invention be limited only by 30 the scope of the appended claims.